

Healthcare  
ent & Unique Healthcare

Safe, Efficient & Unique Healthcare



New Products



Safe, Efficient & Unique Healthcare

Head off: 10 Abo Sholoaa  
Sehet Elzاهر, Cairo, Egypt  
Tel. : ( +202 ) 2593 14 59  
Fax : ( +202 ) 2593 14 67  
[www.abpharma.com.eg](http://www.abpharma.com.eg)

# Purell®


Stay Well With Purell



- Kills 99.99% of most common germs that may cause illness.
- Works in as little as 15 seconds.
- No water or towels needed.
- Specially formulated with 62% ethyl alcohol, moisturizers and vitamin E to leave hands feeling sanitized, soft and refreshed.

Alcohol Hand Gel  
Kills 99.99 %  
of microorganisms

America's No. **1** instant hand sanitizer



Meets FDA's TFM for hand sanitizers





# HMTS

LOW TEMPERATURE PLASMA STERILIZER

Looking for Excellence

- Shortest sterilization time
- Lowest operating cost
- Strongest penetration
- Lowest temperature
- Safest sterilizer



## Accessories



Tray



Incubator



Sterilant Agent



Pouch Roll




C.I. Tape



C.I. Strip



SCBI



Non-Woven Sheet

## The Best Choice



SES  
HMTS-80  
HMTS-120

### **S** Strongest penetration

HMTS is best suited for the sterilization of micro lumen.  
HMTS has expanded the range of sterilization.  
HMTS has powerful penetrability, which increases the sterilizing.  
HMTS sterilization efficacy is not affected by the plasma sheath.

### **L** Lowest operating cost

HMTS has lower installation & operating expenses.  
HMTS reduced management & maintenance cost.

### **S** Safest low temperature sterilization

HMTS final by-products are only oxygen and water vapor.  
HMTS complies with all international environment laws and regulations.

### **S** Shortest sterilization time

HMTS has fast operation time 35 - 40 min.

### **L** Lowest temperature

35° C - 40° C  
HMTS avoids damage & discoloration to medical equipments  
HMTS extends expensive medical equipment's functional life.

### **E**asiest to use

HMTS has touch screen LCD monitor.  
HMTS is available in arabic software  
HMTS installation is simple-plug it in.  
HMTS has printer & 32 MB memory card.







## Gerson N95 2130 Mask

### Properties

- 3D cup shaped respirator folds flat.
- "Away-from-the-face" Design for comfort & non "muffled" conversation.
- Unique side "panels" provide added structural durability.
- Double way protection (protects both the patient and health care worker)
- Synthetic blood penetration resistance as per ASTM F 1862-00a
- Low Breathing resistance as per NIOSH 30 CFR 11 section 11.140-9.
- No skin sensitivity properties. Code A20125f, PSRR # 2002-096.
- FDA Approved , 510 (K) Number: K050193, Regulation Number : 878.4040
- fit well , Equivalent to CE model.  
(Nose clip, elastic band and well sealed edges)



## Gerson Smart\*Mask

### Properties

- 3D cup-shaped respirator with
- "Away-from-the-face" Design for comfort & non "muffled" conversation.
- Ultra-Soft Inner Layer.
- Soft Closed-Cell Foam Nosepiece.
- Durable, Ultrasonically Welded Non-Latex Head Bands.
- No Uncomfortable Pressure Points.
- Adjustable Integrated Nosepiece.
- Flexible Contoured Edges
- Contribute to Better Fit.
- Small Size Available
- NIOSH approved, 42 CFR 84



**NEW!**  
**Comfort  
MAXX**  
Super Valve



40% lower breathing resistance than traditional respirators



# Protect Yourself & Your Patients



## Face Mask with Shield

### BE CAREFUL

the risks of eye contamination are often not sufficiently recognized

### Properties

- 99% bacterial filtration efficiency.
- Shield sonic welded to mask.
- Anti-fog face shield.
- Easy breathing.
- Comfortable with or without eye wear.
- Fluid resistant mask.
- Tie-on mask with shield also available.



## Autoclavable Clogs

### Properties

- Autoclavable.
- Anti-microbial.
- Anti-static.
- Healthy.
- Hygienic.
- Comfortable.
- Superior quality.





# Tristel

Chlorine Dioxide

*Sporocidal in few seconds*

*New Dimensions of high level disinfection*



Surface



- Work surfaces
- Beds & mattress covers
- Instrument tables
- Monitor cases
- Control panels
- Key boards

Ultrasound



Instruments



- Rapid Action Time 5 min.
- Safe To User
- Economic
- Error-free System
- Single Use
- Corrosion Inhibitors
- Biodegradable 100%

Surface



The world's 1st rapid action sporocidal foam

Tristel is the most widely used sterilising solution in UK Hospitals



# Keep Your System Standard VALIDATION

## For Accreditation



A.B.Pharma has the qualifications, equipments & well trained staff to provide standard validation of clean areas in hospitals, pharmaceuticals, food & beverages according to all international standards. Our Validation licensed from CZECH Co. which is permanent member of the European standardization committee & provider of clean room validations according to standards from ISO CE GMP TUV FDA

## System Sure II

### Hygiene Monitoring System

#### Validation Tests

- HEPA Filter Leak Test
- Air Particles Counting
- Air Flow Differential Pressure
- Air Flow & Velocity Measurements
- Air Change Per Hour (ACH)
- Air Flow Visualization Display Analysis
- Temperature & Humidity
- Noise & Light Lux Measurements
- Recovery Test
- Microbial analysis



#### Applications

Validation of Air Quality System in :

- Hospitals  
Operating Theaters, ICUs & Labs
- Pharmaceuticals and Clinical Pharmacy  
Labs, Clean areas & Sterile Areas
- Food & Beverages  
Labs
- Laboratories & Research Centers  
Biological Safety Cabinets & Hoods



#### Where Rapid Detection Of Contamination Is Crucial

- Small, lightweight
- Easy to use, extremely sensitive and affordable
- Detects contamination rapidly: 1.5 Sec.
- Enables immediate corrective action and avoid recalls
- Provides accurate and traceable verification of the hygienic status of surface & water.
- Validates HACCP and SSOP programs.

#### Datasure II software

- Download results to a PC
- Produce color trend analysis graphs & reports
- Identify problem areas and manage HACCP & SSOP tracking requirements





L a m i n a r A i r F l o w

Safe Surgery Saves Lives

O P i

Island Of Ultra Clean Air

- 4m<sup>2</sup> Sterile Area.
- Decreases SSI By 52%.
- Highest Standard Of Air Quality, ISO 5.
- Maintains Sterility Of Surgical Instruments.
- Decreases The Settlement Of Microbes On The Surface By 75%.
- Reduces Patient's Stay By 2.5 Times.
- Heavy Worker, 24 Hours.
- Maintenance Free.
- Self Control Computerized System.

Renovate Your Old Operating Room To Meet ISO 5 Standard